

NOR-LAKE[®] SCIENTIFIC

DATA LOGGER / THERMOMETER / MONITOR WITH DIGITAL DISPLAY WITH TEMPERATURE ALARM

THERMALARM WITH USER-DEFINED TIME INTERVAL

- Alarm tracking stores in memory when each alarm event occurred and duration
- Meets current CDC requirements for vaccine storage and monitoring
- Rolling memory structure, maintains most recent 525,600 temperature observations
- Stores up to 10 unique alarm events
- Memory can be cleared on unit after downloading data to USB stick
- Hassle-free retrieval of data - unit can remain in use while downloading and analyzing data
- Status indicators—low battery, memory full, USB data transfer and active alarm state
- Supplied: stand, Velcro[®], magnetic strips, wall mount, Traceable[®] Certificate, batteries



Model	Range	Resolution	Accuracy	Probe
ThermAlarm	-50.00 to +70.00°C	0.01°	± 0.25°C	1 bottle

Nor-Lake[®] Scientific

727 Second Street • Hudson, WI 54016
800-477-5253 • www.norlakescientific.com

