

NOR-LAKE® SCIENTIFIC ENVIRONMENTAL ROOMS

ENVIRO-LINE WALK-IN ROOMS from Nor-Lake Scientific reproduce and closely monitor environments with a precise combination of humidity and temperature. Choose a multi-temperature or constant temperature room. Environmental Rooms are customized to meet your testing and storage needs.

- Suitable for Research, Biological Studies, Shelf Life, Stability Testing, Dry Rooms, Blood & Plasma Storage, General Life Science Applications and Working Laboratories
- Custom Rooms Available to Suit Specific Storage and Application Requirements
- Modular Construction for Easy Installation and Relocation
- Available Microprocessor Controllers, Alarms, Monitors, Temperature and Humidity Chart Recorders, Shelving, Ramps, Lighting, Humidification and Dehumidification



MINI-ROOM WALK-IN CONTROLLED TEMPERATURE CHAMBERS

- Single Room & Combination Cold Rooms, Warm Rooms and Freezers
- Temperatures +4°C, -20°C, -30°C, +2°C to +55°C
- 10,000 Models to Select From
- Installation is Fast and Easy
- Wall or Ceiling Mount Capsule Pak™ Refrigeration System Keeps all Components Outside Allowing More Storage Inside
- Standard Sizes from 4' x 5' to 8' x 14'
- Indoor and Outdoor Models
- Options: Microprocessor Control, Chart Recorder, Humidification, Dehumidification, Duplex Outlets, Shelves



NORLAKE SCIENTIFIC



 Committed to our environment

NOR-LAKE® SCIENTIFIC

727 Second Street
Hudson, WI 54016

800-477-5253

715-386-2323

715-386-4290 FAX

www.norlakescientific.com

NOR-LAKE[®] SCIENTIFIC

LABORATORY AND SCIENTIFIC EQUIPMENT

For Life Science, Hospital, Clinical, Biomedical, Scientific and Industrial Applications

Nor-Lake[®] Scientific Select[™] Ultra-Low Freezers

- Upright Models Available in 14.3, 18.9, 23.3 & 27.5 ft³ Internal Storage Capacity
- Chest Models Available in 4.6 ft³ Internal Storage Capacity
- Operating Temperature Range -50°C to -86°C
- Next Generation 1HP Low-Noise, High-Efficiency Compressors
- Exclusive Cascade Refrigeration and Cabinet Design
- Programmable Logic Microprocessor Controller with LCD Digital Display
- High/Low Audible and Visual Temperature Alarms, Remote Alarm Contacts
- Door Ajar Alarm, On Board Diagnostics, Surge Protector, Password Protection, Real Time Clock, Battery Backup of Controller Display
- Available Options: CO₂ and LN₂ Backup Systems, Chart Recorder, Inventory Storage Systems
- UL, C-UL Listed



Nor-Lake[®] Scientific Select[™] and Premier[™] Laboratory & Pharmacy Refrigerators and Freezers

- Available in 1, 2 and 3 Glass or Solid Door Models with 5, 24, 33, 52 & 80 ft³ Storage Capacity
- Pass-Thru: Chromatography and Sliding Door Models Available
- Microprocessor Controller, Digital Display, Hi/Lo Audible & Visual Alarms with Remote Alarm Contacts, Key Lock Doors, Shelving, Casters
- Options: Internal Electrical Outlet, Temperature Recorder, Stainless Steel Interior and Exterior, Access Ports, Drawers
- UL, C-UL Listed



Nor-Lake[®] Scientific Blood Bank Refrigerators Plasma Freezers

- Meets AABB, ANRC and FDA Blood Storage Requirements
- Available in 1, 2 and 3 Glass Door Models with 24, 33, 52 & 80 ft³ storage capacity
- Microprocessor Controller, LCD Digital Display (Air and Product) Password Protection, Hi/Lo Audible & Visual Alarms with Remote Alarm Contacts, Temperature Chart Recorder, Drawers, Glass Doors, Key Lock Doors, Casters
- UL, C-UL Listed



Nor-Lake[®] Scientific Stability Test Chambers BOD Refrigerated Incubator CO₂ Incubators Lighted Chambers Warming Cabinets

- Various Models with Temperature only or Temperature and Humidity Control
- Available with 1, 2, and 3 Doors with 24, 33, 52 & 80 ft³ storage capacity
- Programmable Microprocessor Controller, LCD Digital Display, Hi/Lo Audible & Visual Alarms with Remote Alarm Contacts, Key door lock, casters, Shelves, Stainless Steel Interior
- Options: Chart Recorder, Access Port, Internal Electrical Outlet, Extra Shelves, RS485, Stainless Steel Exterior
- UL, C-UL Listed



Nor-Lake[®] Scientific Table Top and Undercounter Refrigerators, Freezers, Combination Refrigerator/Freezer

- Models Available for +4°C, -20°C & -30°C Temperature Ranges for General Purpose, Blood, Plasma, Reagent and Vaccine Storage
- Auto defrost and manual defrost
- Microprocessor controls with alarm



NOR-LAKE[®] SCIENTIFIC

727 Second Street
Hudson, WI 54016

800-477-5253

715-386-2323

715-386-4290 FAX

www.norlakescientific.com

Part Number: 095406

Revision Date: 05/18

