


NOR-LAKE[®] SCIENTIFIC

GENERAL PURPOSE LABORATORY

REFRIGERATOR/FREEZER UPRIGHT AUTO DEFROST

Standard Features:

- Operating temperature: +4°C refrigerator, -10°C freezer
- Auto defrost, forced air
- Adjustable temperature control, independent for refrigerator and freezer
- 3 adjustable sliding full width shelves (refrigerator)
- 1 adjustable shelf (freezer compartment)
- 3 independent removable clear plastic drawer storage bins (refrigerator)
- 2 full-width door storage bins (freezer)
- 3 full-width fixed door bins (refrigerator)
- 1 door storage bin with removable hinged clear plastic cover (refrigerator)
- Interior light (refrigerator)
- 2 adjustable leg levelers
- Automatic door closer
- White interior, exterior, door and handle
- Complete cabinet includes CFC-free high density polyurethane foam insulation
- Warranties: 12 months parts and labor (US and Canada), 12 months parts (International)



Digital Thermometer Alarm (optional)

Optional Features:

- Data Logger/Thermometer Monitor with Digital Display with alarm (field install)
- Export crating



LRF201WWW/0



LRF201WWW/0


NORLAKE[®] SCIENTIFIC




Committed to our environment

NOR-LAKE SCIENTIFIC
A Division of
NOR-LAKE, INCORPORATED
727 Second Street


Hudson, Wisconsin 54016
800-477-5253
715-386-2323
800-388-5253 PARTS
715-386-4290 FAX
www.norlakescientific.com




TOP VIEW



INNER DOOR VIEW



FRONT AND SIDE VIEWS



SHELF VIEW

SPECIFICATIONS	
Model	LRF201WWW/0
Crated Weight (lbs) (kg)	240 (108)
Crated Height (in) (cm)	74 (188)
Crated Width (in) (cm)	42 (106)
Crated Depth (in) (cm)	42 (106)
Interior Height (in) (cm)	
Refrigerator Compartment	39-1/2 (100)
Freezer Compartment	17-1/4 (43)
Interior Width (in) (cm)	
Refrigerator Compartment	26 (66)
Freezer Compartment	24 (60)
Interior Depth (in) (cm)	
Refrigerator Compartment	26-1/2 (67)
Freezer Compartment	22 (55)
Overall Height (in) (cm)	69 (175)
Overall Width (in) (cm)	30 (76)
Overall Depth (in) (cm)	34 (86)
Gross Cubage (CuFt) (L)	
Refrigerator Compartment	15.0 (425)
Freezer Compartment	5.5 (156)
Number of Shelves	4
Shelf Area (SqFt) (Sqm)	12.1 (1.1)
Refrigerant	R-134a
Total Amp Draw	6

Specifications subject to change without notice.

Voltage Model Suffix Code	Voltage Description	NEMA Plug	NEMA Receptacle
/0	115/1PH, 60HZ	5-15P	5-15R



NORLAKE SCIENTIFIC

Revision Date: 08/15
©2015 NOR-LAKE, INC.
Printed in the U.S.A.
Part Number: 139507

